

FOAMING CLEANER

Clinging Foam Disinfectant

Form: Aerosol Liquid

Color: Colorless

Odor: Lemon

pH: 12-13

Solubility (in water): Soluble
(in mineral spirits): No data available

VOC Content (% by weight): 13

Flash Point: None

Specific Gravity: 1.00 g/cm³

Density: 8.35 lbs/gal

Storage Stability (at 70°F): 1 year

<u>Ingredients</u>	<u>C.A.S. #</u>
Water	7732-18-5
Sodium Metasilicate	6834-92-0
Diethylene Glycol Monobutyl Ether	112-34-5
Propellant.....	Blend
Alkyl (60%C14, 30%C16, 5%C12, 5%C18)	
dimethylbenzyl ammonium chloride.....	63449-41-2
Alkyl (68%C12, 32%C14) dimethylethylbenzyl	
ammonium chloride	85409-23-0

Ideal for...
• Hospitals
• Clinics
• Restrooms
• Schools
• Hotels

DIRECTIONS:

Disinfection, Deodorization and Cleaning: Shake can before using. Hold container six to eight inches from surface to be treated. Spray area until it is covered with white foam. Allow foam to penetrate and remain wet for 10 minutes. No scrubbing is necessary. Wipe off with a clean cloth, mop or sponge. This product will not leave grit or soap scum.

For heavily soiled areas, a pre-cleaning step is required before treatment.

SEE PRODUCT LABEL FOR DETAILED DIRECTIONS.

HMIS®			NFPA®		
	Severe	4	Extreme		
Health.....	Serious	3	High	Health.....	1
Flammability	Moderate	2	Moderate	Flammability	4
Reactivity	Slight	1	Slight	Reactivity	0
Personal Protection	Minimal	0	Insignificant	Special Precautions	None